



# DAMSystem Product Price List



DAM2

## DAM2 Drosophila Activity Monitor


- 32 tubes, 5 or 7mm diameter
- Dual consolidated IR beams per tube
- Available models
  - DAM2 for 5mm tube = \$485.
  - DAM2-7 for 7mm tube = \$545.
- DAM2-HS for high speed counting of flies pulled by vacuum = \$600. (6 channels only, 5mm tubing)



DAM5M

## DAM5M Drosophila Activity Monitor


- 32 tubes, 5 mm diameter
- 4 independent IR beams per tube
- Reports Counts, Moves, Dwell, Position, Rest
- Built-in clips hold tubes securely
- Available models
  - DAM5M = \$785.



LAM25H

## LAM10/16/25 Locomotor Activity Monitor

- 32 tubes, 10, 16, or 25mm diameter
- 3 or 9 (high res) consolidated IR beams per tube
- Available models
  - LAM10, LAM16, LAM25 = \$875.
  - LAM10H, LAM16H, LAM25H = \$1,175. (high res)
  - LAM10H-3, LAM16H-3, LAM25H-3 = \$3200. 3-board stack monitors 3 axial positions per tube
- LAMVTS Vertical Tube Support plate = \$75.



LAM50H

## LAM50H Locomotor Activity Monitor


- 32 tubes, 50mm diameter
- 9 consolidated IR beams per tube
- Available models
  - LAM50H = \$1,350.

**TRIKINETICS**



### DPM Drosophila Population Monitor

- Single 25mm vial for group activity monitoring
- 3 IR beam sets monitor activity at 3 separate axial locations
- 15 beams per set robustly detect individual *d. melanogaster* walking on the tube walls
- Operates with horizontal or vertical vial
- Available models
  - DPM = \$750. (includes 25x95mm glass shell vial)



### DEM Drosophila Eclosion Monitor

- Counts falling flies after eclosion from pupae disk
- Solenoid tapper periodically dislodges new imagoes from disk and funnel surfaces
- Available models
  - DEM = \$1,750. (includes funnel and pupae disk)
- Accessories:
  - FUN100 100mm funnel = \$30.
  - PDSK Pupae disk = \$18.



### MB5 MultiBeam Activity Monitor

- 16 tubes, 5, 7, or 10mm diameter
- 17 independent IR beams per tube report Counts, Moves, Dwell, Position, Rest at up to 1hz
- Position preference vs food, odor, light, etc
- Available models
  - MB5, MB7, MB10 = \$1,550.
- PSIU24 Power Supply Interface Unit = \$350. supports 32 MB monitors, includes all cables



### DEnM Drosophila Environment Monitor

- Measures incubator conditions alongside activity data for archival record
- Precision measurements verify incubator setpoints
- Temperature: 0-70C, +/- 0.1 degC
- Relative Humidity: 5-95%, +/- 2%
- Ambient Light: 0-2500 lux, +/- 5%
- Available models
  - DEnM = \$450. (Temp, RH, Light)
  - DEnM-P = \$450. (Barometric Pressure, in Hg)




### PSIU9 Power Supply Interface Unit

- Connects data collection computer to activity monitors. One required per system.
- Includes PS9-1 AC power supply, AC line cord, USB cable, and 10' extension data cable with coupler
- Available models
  - PSIU9 = \$350. (includes PS9-1)
  - PS9-1 Power supply = \$30. (9VDC, 1A)



### LC4 Light Controller

- On/off switching of 4 AC outlets
- Pulse duration from 0.1 second to 99 hours
- Single, periodic, or random pulse repetition
- Incubator lighting control
- Vortexer pulsing for sleep deprivation
- 3A, 240VAC max for each outlet, 50/60 hz
- Available models
  - LC4 = \$750.



### MAN2 Gas Distribution Manifold


- Facilitates gas flow thru DAM2 5mm tubes
- Tubes snap into manifold via gasket seal
- Exit holes in plastic tubes (optional)
- Disassembly for cleaning
- Available models
  - MAN2 for 5mm tube = \$400.
  - PPT5x65D2E (50) = \$50. drilled PC tubes, 2x.030" holes, 1cm to end



### DFM Drosophila Funnel Monitor

- Counts falling flies under inibriometer column
- 75mm funnel diameter
- 8mm OD funnel neck
- Available models
  - DFM = \$800. (includes funnel)





Vortexer and VMP

### Vortexer Mounting Plate


- Fits Troemner 945007/945008 vortexer
- Securely holds activity monitors during shake cycle
- Compatible with LC4 Light Controller for precision computerized sleep deprivation
- Includes all necessary hardware

### Mounting Plates


- VMP = \$575. (holds 4 DAM2 monitors)
- VMP-UD = \$250. (upper deck for VMP)
- VMP-LAM = \$575. (holds LAM10/16/25)
- VMP-MB = \$750. (holds 2 MB5)

### Troemner vortexer with VMP Mounting Plate installed


- TVOR-120 = \$3,000. (120V = 945007)
- TVOR-230 = \$3,000. (230V = 945008)




VMP and VMP-UD holds 8 DAM2 units



VMP-LAM




VMP-MB



### Monitor Tubes (outside diameter x length, mm)

- PGT5x65 Pyrex Glass (100) = \$70.
- PPT5x65 Polycarbonate (100) = \$15.  
PPT5x65 Polycarbonate (1000) = \$120.  
PPT5x65 Polycarbonate (7500) = \$675.
- PGT5x80 Pyrex Glass (100) = \$75.
- PPT5x80 Polycarbonate (100) = \$16.
- PGT7x65 Pyrex Glass (100) = \$80.
- PGT10x100 Pyrex Glass (100) = \$100.
- PGT16x100 Pyrex Glass (43) = \$50.  
PGT16x100 Pyrex Glass (250) = \$230.
- PGT25x125 Pyrex Glass (93) = \$180.
- PGT60x300 Pyrex Glass (1) = \$50.


Tubes are open on both ends; glass are fire glazed.



### Tube Caps

- CAP5-Black (100) = \$10. (fits 5mm tube)
- CAP5-Yellow (100) = \$15.
- CAP7-Black (100) = \$11.
- CAP10-Black (100) = \$12.
- CAP16-Black (100) = \$15.
- CAP25-Black (100) = \$20.

Black caps are vinyl plastic, yellow are PVC plastic. Both may be cut to length with scissors and drilled for ventilation.



### Cables and Connectors

- SPLT5 5-way splitter = \$10.
- CPLR Coupler = \$5.
- CAB2 2' monitor cable = \$5.
- CAB10 10' extension cable = \$8.
- SPLT8-6 8-way splitter for multibeam = \$80.
- CAB2-6 2' cable, 6 cond for MB = \$6.
- CAB10-6 10' cable, 6 cond for MB = \$9.

# Ordering Notes

How to order	Send a Purchase Order by email to sales@trikinetics.com, including item list, shipping address, and email address for invoice submission. Credit card payment is accepted with an added 3% (domestic) or 4.5% (foreign) PayPal handling fee. A PayPal account is not required.
Basic System	A basic system will include the PSIU9, one or more activity monitors, and a customer-supplied Macintosh or Windows PC with USB port for data collection and storage.
Tubes and Caps	Unless otherwise noted, monitors do not include tubes or caps, which may be ordered separately.
Wiring Accessories	All systems are shipped with appropriate DAMSystem network wiring accessories, including cables, couplers, and 5-way splitters. Monitor cables are 50cm (2' ) long and extension cables are 3m (10' ) long unless otherwise specified.
Monitor Numbers	All units are identified by a unique numeric address in the range 1:120 which must be specified at the time of order. Addresses should not be duplicated within a system, and will be assigned consecutively from #1 unless otherwise specified.
Software	The DAMSystem data acquisition software package and PSIU9 USB driver are available for download at no charge from the TriKinetics web site, www.trikinetics.com.
Modifications	Custom variants of any and all of the above-listed products and accessories are available on request.
Foreign Orders	Overseas shipment is available for any item. Foreign customers should be advised that VAT or other import duty will likely be assessed on arrival, and may slow delivery if preparations are not made in advance.
Suggestions	Comments on product performance, and suggestions for product improvement, are particularly important to us, and always welcome.
Warranty and Repair	Defective units will be repaired at no charge.
Specifications	All specifications are nominal, and subject to change without notice to improve product reliability or performance. No liability is assumed for loss of data or consequential damages, and all trademarks are the property of their respective owners.
Acknowledgments	These units were developed with the continuing assistance and support of scientists in many laboratories around the world, and this support is gratefully acknowledged. We are especially indebted to Prof. Michael Rosbash and Prof. Jeffrey Hall at Brandeis University for their cooperation and assistance in the early days.
Contact	TriKinetics Inc Waltham, MA USA sales@trikinetics.com ©1986-2018